



affinisep

The Art of making
sample preparation easier



PROTEOMICS
METABOLOMICS
BIOLOGICAL APPLICATIONS



SPE Membrane	5
SPE membrane technology for microelution	6
Focus on SPE membrane microelution tools	8
Our product ranges for proteomics and biomolecular applications	11
BioSPE™	12
BioSPE™ ready-to-use kits for proteomics and metabolomics...	13
BioSPE™ tools – Capacity	14
BioSPE™ PurePep for peptide desalting	17
BioSPE™ PepFrac for peptide fractionation	21
BioSPE™ PureProt for intact protein purification	24
BioSPE™ GlycaClean for glycopeptide enrichment	26
BioSPE™ BeadRem for magnetic beads removal	28
Choose the perfect BioSPE™ product for your proteomic application	33
BioSPE™ PureOmics for metabolite purification	35
Choose the perfect BioSPE™ product for your metabolomic application	38

AttractSPE® Disks	39
AttractSPE® Disks chemistries	40
AttractSPE® Disks tools capacity	41
AttractSPE® Disks Tips – StageTips 200µL and 1mL	42
AttractSPE® Disks Tips – StageTips 10µL	43
AttractSPE® Disks Tips C18: NEW ECO PACKAGING!	44
AttractSPE® Disks 96 well plates for microelution	46
AttractSPE® Disks 96 well plates – 1mL	47
AttractSPE® Disks Spin Columns	50
AttractSPE® Disks Cartridges	53
FAQ about AttractSPE® Disks and BioSPE™	55
Membrane filter & powder-based SPE	60
Removal of proteins & lipids	61
AttractFiltr membrane filter	63
Powder-based SPE cartridges and 96 well plates	64

SPE Membrane

Our latest innovation for
microelution in proteomics and
biomolecular applications



SPE membrane technology for microelution

Our innovative SPE membranes are made of very small sorbent beads tightly embedded, thus leading to **thin, dense, soft and uniform membranes with high capacity** for extraction/separation, purification and concentration of molecules.


Thanks to their unique advantage, our SPE membranes are useful for purification of very small sample volumes in proteomics, genomics, metabolomics, bioanalysis, biomarker discoveries and biological applications. They are applied for spinnable, automatable, high throughput microelution and nano elution.

These SPE membranes offer outstanding sample preparation efficiency and reproducibility of results. Since the diffusion distance between particles is minimized, adsorption is more efficient, and extraction can be accomplished using **very low sorbent mass**.

These properties are giving to our SPE membranes a significant improvement of mass transfer kinetics compared to traditional packed SPE particles. As monoliths, our SPE membranes are self stand and require no frits for immobilizing the column bed (unlike traditional SPE products), allowing 100% recovery of the original sample volume.

Advantages of SPE membranes

- No need to frits or filters
- Reduced dead volume
- Small elution volumes
- High sample recovery
- Reduced time for eluate evaporation
- Higher throughput
- Channeling effects eliminated
- Excellent reproducibility
- Concentration of the sample



SPE membrane technology for microelution

Our SPE membranes can be used to miniaturize SPE for small operating volume of fluid samples. **4 formats** are **available** for microelution. Thanks to the use of the SPE membranes for all formats, **the change of format or the scale up of the process is easy.**

Capacity of molecules of interest increases

□ StageTips (SPE tips) and 96 StageTips well plates for microelution

μ SPE columns designed by immobilizing a uniform membrane inside a pipette tip (StageTips per unit or as 96 well plates)



□ 96 SPE Well Plates

SPE 96 well plates designed by immobilizing a uniform membrane



□ SPE Spin Columns

SPE columns designed by immobilizing a uniform membrane inside a centrifuge SPE tube



□ SPE Cartridges

SPE cartridges designed by immobilizing a uniform membrane



Concentration of molecules of interest increases


StageTips or SPE Tips, SPE Spin columns, 96 SPE well plates and SPE cartridges are tools for proteomic, genomic and bioanalysis (small molecules cleanup). In order to give the most exhaustive applications, they are proposed with a large variety of formats and capacities.



**SPE Tips
StageTips**



**SPE Spin
columns**



**96 SPE well plates for
microelution**



SPE Cartridges



Advantages of SPE membrane tools

- **Remove interfering contaminants** — significantly reduce signal suppression and improve signal-to-noise ratios and sequence coverage
- **Simplify optimization** — processing yields high-quality spectra and is effective for a variety of reverse-phase-compatible contaminants
- **Robust** — work with a wide variety of load volumes and concentrations; no need to reduce sample volume before application
- **Convenient** — easy to handle and require no special equipment to process multiple samples simultaneously (unlike tip-driven systems that require one sample to be processed at a time)
- **Sensitive** — special resins allow excellent recovery percentages, even at low (sub-picomole) sample loads

Capacity

Each format is available with **different capacities**, that is to say different layer thicknesses of the SPE membrane immobilized inside the microelution tool, to better adapt to your needs. For a more reproducible product, each tool contains only **ONE** layer of SPE membrane.





Illustration of the different capacities available for our SPE Tips



As manufacturers of our own SPE membranes, we can offer **custom-made solutions** for complex applications (e.g. stacking of several sorbents with different capacities). Contact us for more information.

Use of SPE membrane microelution tools

SPE Tips or StageTips, SPE Spin columns and 96 SPE well plates are easily spinnable tools.



96 SPE well plates and **SPE cartridges** can be used with positive pressure and/or vacuum manifolds and are compatible with SPE automates.



DID YOU KNOW?

96 SPE well plates can also be used by centrifugation. Partial use is possible, without consequences for the unused wells.

Elution of samples from our **SPE Tips** can be performed by **centrifugation** (with adapted devices) or with the use of a positive pressure or vacuum manifold. **SPE Tips** can also be used with automates, with a refilling by the top of the tips.

Our **SPE Tips** were not developed for a use by pipetting. Thus, we strongly recommend to use our **SPE Tips** by centrifugation.

Our product ranges for proteomics and biomolecular applications

Affinisep offers **two ranges of products** for proteomics and biomolecular applications, **both based on our innovative SPE membrane technology**.

BioSPE™

A new range of **ready-to-use kits** specially developed and designed for **various applications in bottom-up proteomics** (peptide desalting, peptide fractionation, removal of magnetic beads after SP3 procedure, enrichment of glycans or glycopeptides), **top-down proteomics** (purification of intact proteins) or **metabolomics** (purification of small molecules).



BioSPE™ is the perfect product range if you are looking for **turnkey SPE solutions for well-defined applications!**


AttractSPE®Disks

A comprehensive product range based on sorbent chemistries, with a **wide variety of silica-based and polymeric sorbents**, for extraction and purification of analytes in biological samples such as urine, plasma, serum, saliva....



AttractSPE®Disks is the perfect product range if you are looking for **SPE solutions with specific sorbent chemistries** or for **complex applications!**

BioSPE™



Ready-to-use kits
for

Peptide desalting
Peptide fractionation
Intact protein purification
Glycopeptide enrichment
Magnetic beads removal
Metabolite purification


BioSPE™ product range consists of **ready-to-use kits** for the most common applications in proteomics and metabolomics workflows.

Supplied as **turnkey solutions**, these kits are composed of:

- ❑ SPE columns for a given application such as desalting, enrichment, fractionation...
- ❑ Instructions of use and detailed protocol specific to a given application
- ❑ Adapters and/or collection tools (each kit contains collection tools for the conditioning/loading/washing fractions and for the elution fraction)
- ❑ Some reagents

BioSPE™ 96 well plates are supplied with **1mL collection plates** and **BioSPE™ Spin columns** are supplied with **2mL collection tubes**.

BioSPE™ Tips are supplied with **spin adapters and collection tubes** for their centrifugation. Very easy to use, these adapters are compatible with any 1.5 and 2mL centrifugation tube (in the case of 10µL and 200µL tips) and 15mL centrifugation tube (in the case of 1mL tips).



Products from the BioSPE™ range are available with **two different capacities, which are referred S (Standard) and H (High)** in the product designations. This designation refers to the layer thickness of SPE membrane. For a more reproducible product, each tool contains only **ONE** layer.








Illustration of the two capacities available for our BioSPE™ Tips

Following references (p17 to p27) will have **XX** indications, which have to be replaced by the capacity you require (S or H).

To help you select the formats that are the most adapted to your applications based on your sample amounts, the following table lists the **loading capacities** (maximal peptide/protein/metabolite amounts that can be retained) of our different **BioSPE™** formats available.

Products		Capacity (µg)	
		Standard (S)	High (H)
	BioSPE™ Tips - 10µL	7	25
	BioSPE™ Tips - 200µL	15	50
	BioSPE™ Tips - 1mL	35	105
	BioSPE™ 96 well plate for microelution	35	105
	BioSPE™ 96 well plate 1mL	500	1500
	BioSPE™ Spin - Micro	200	600
	BioSPE™ Spin - Mini	300	900
	BioSPE™ Spin - 15mL	1500	4500
	BioSPE™ Spin - 50mL	3000	9000
	BioSPE™ Cartridges – 1mL	250	750
	BioSPE™ Cartridges – 3mL	700	2100
	BioSPE™ Cartridges – 6mL	1500	4500



Application: Clean-up of **peptide mixtures** after **enzymatic digestion** in bottom-up proteomics to **remove salts** that can interfere with peptide ionization/detection and damage mass spectrometers



Advantages

- **Up to 97% proteins identified** after desalting of HeLa digest on **BioSPE™ PurePep**
- **Good recovery of peptides on the whole polarity range**
- **Excellent repeatability (RSD < 3%)**



DID YOU KNOW?

For desalting of the **most hydrophilic peptides**, we recommend the use of **BioSPE™ PurePep Broad**.



Advantages

- **Efficient retention of small polar peptides**
- **Stable under extreme pH conditions (pH 0-14)**
- **BioSPE™ PurePep Broad can run dry** during SPE process

BioSPE™ PurePep for peptide desalting

Designation	Kit composition	Reference*
BioSPE™ PurePep – Tips 10µL Single Cell & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PurePep.SC.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PurePep.SC.10.96
BioSPE™ PurePep – Tips 10µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PurePep.XX.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PurePep.XX.10.96
BioSPE™ PurePep – Tips 200µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PurePep.XX.200.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PurePep.XX.200.96
BioSPE™ PurePep – Tips 1mL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PurePep.XX.1000.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PurePep.XX.1000.96
BioSPE™ PurePep – Micro Spin	Spin columns (50/pk) + collection tubes (100/pk)	µSpin-PurePep.XX.50
BioSPE™ PurePep – Mini Spin	Spin columns (50/pk) + collection tubes (100/pk)	Spin-PurePep.XX.50
BioSPE™ PurePep – Spin 15mL	Spin columns (25/pk) + collection tubes (50/pk)	Spin15-PurePep.XX.25
BioSPE™ PurePep – Spin 50mL	Spin columns (25/pk) + collection tubes (50/pk)	Spin50-PurePep.XX.25

*Replace XX by the capacity (S or H, see p15) of your choice. e.g.: *KT-Adapt-Tips-PurePep.S.200.96*

Designation	Kit composition	Reference*
BioSPE™ PurePep – 96 plate for microelution & collection plates	Microelution 96 SPE Well Plate (1/pk) + collection plates (2/pk)	KT-Coll-μ96W-PurePep.XX.1
BioSPE™ PurePep – 96 plate 1mL & collection plates	1mL 96 SPE Well Plate (1/pk) + collection plates (2/pk)	KT-Coll-96W-PurePep.XX.1
BioSPE™ PurePep – Cartridges 1mL	1mL SPE cartridges (50/pk)	CAR1-PurePep.XX.50
BioSPE™ PurePep – Cartridges 3mL	3mL SPE cartridges (50/pk)	CAR3-PurePep.XX.50
BioSPE™ PurePep – Cartridges 6mL	6mL SPE cartridges (50/pk)	CAR6-PurePep.XX.50

*Replace XX by the capacity (S or H, see p15) of your choice. e.g.: CAR1-PurePep.S.50

BioSPE™ PurePep Broad for peptide desalting on broad pH range

Designation	Kit composition	Reference*
BioSPE™ PurePep Broad – Tips 10µL Single Cell & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PurePep- Broad.SC.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PurePep- Broad.SC.10.96
BioSPE™ PurePep Broad – Tips 10µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PurePep- Broad.XX.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PurePep- Broad.XX.10.96
BioSPE™ PurePep Broad – Tips 200µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PurePep- Broad.XX.200.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PurePep- Broad.XX.200.96
BioSPE™ PurePep Broad – Tips 1mL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PurePep- Broad.XX.1000.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PurePep- Broad.XX.1000.96
BioSPE™ PurePep Broad – µSpin	Spin columns (50/pk) + collection tubes (100/pk)	µSpin-PurePep- Broad.XX.50
BioSPE™ PurePep Broad – Spin	Spin columns (50/pk) + collection tubes (100/pk)	Spin-PurePep- Broad.XX.50
BioSPE™ PurePep Broad – Spin 15mL	Spin columns (25/pk) + collection tubes (50/pk)	Spin15-PurePep- Broad.XX.25
BioSPE™ PurePep Broad – Spin 50mL	Spin columns (25/pk) + collection tubes (50/pk)	Spin50-PurePep- Broad.XX.25

*Replace XX by the capacity (S or H, see p15) of your choice. e.g.: Spin-PurePep-Broad.S.50

Designation	Kit composition	Reference*
BioSPE™ PurePep Broad – Tips 10µL Single Cell & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PurePep-Broad.SC.10.25
BioSPE™ PurePep Broad – 96 plate for microelution & collection plates	Microelution 96 SPE Well Plate (1/pk) + collection plates (2/pk)	KT-Coll-µ96W-PurePep-Broad.XX.1
BioSPE™ PurePep Broad – 96 plate 1mL & collection plates	1mL 96 SPE Well Plate (1/pk) + collection plates (2/pk)	KT-Coll-96W-PurePep-Broad.XX.1
BioSPE™ PurePep Broad – Cartridges 1mL	1mL SPE cartridges (50/pk)	CAR1-PurePep-Broad.XX.50
BioSPE™ PurePep Broad – Cartridges 3mL	3mL SPE cartridges (50/pk)	CAR3-PurePep-Broad.XX.50
BioSPE™ PurePep Broad – Cartridges 6mL	6mL SPE cartridges (50/pk)	CAR6-PurePep-Broad.XX.50

*Replace XX by the capacity (S or H, see p15) of your choice. e.g.: CAR3-PurePep-Broad.H.50




Application: Reduction of sample complexity with peptide fractionation at basic pH for deep proteome sequencing and quantitative analysis



Advantages

- Increased number of identified proteins
- Minimized fraction overlapping and good distribution of peptides in all fractions
- No storage constraints (dry at room temperature) with long shelf life
- Real flexibility of format and capacity to adapt to all samples



BioSPE™ PepFrac for peptide fractionation

BioSPE™ PepFrac ready-to-use kits are supplied with **0.1% triethylamine (TEA) solution**. Up to 8 fractions can be performed thus the kits are supplied with **nine collection tools** (one for conditioning/loading/washing fractions and one for each of the eight fractions).

Designation	Kit composition	Reference*
BioSPE™ PepFrac – Tips 10µL Single Cell & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (225/pk) + 0.1% TEA (25mL)	KT-Adapt-Tips-PepFrac.SC8.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (864/pk) + 0.1% TEA (25mL)	KT-Adapt-Tips-PepFrac.SC8.10.96
BioSPE™ PepFrac – Tips 10µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (225/pk) + 0.1% TEA (25mL)	KT-Adapt-Tips-PepFrac.XX8.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (864/pk) + 0.1% TEA (25mL)	KT-Adapt-Tips-PepFrac.XX8.10.96
BioSPE™ PepFrac – Tips 200µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (225/pk) + 0.1% TEA (100mL)	KT-Adapt-Tips-PepFrac.XX8.200.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (864/pk) + 0.1% TEA (200mL)	KT-Adapt-Tips-PepFrac.XX8.200.96
BioSPE™ PepFrac – Tips 1mL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (225/pk) + 0.1% TEA (100mL)	KT-Adapt-Tips-PepFrac.XX8.1000.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (864/pk) + 0.1% TEA (200mL)	KT-Adapt-Tips-PepFrac.XX8.1000.96
BioSPE™ PepFrac – Micro Spin	Spin columns (50/pk) + collection tubes (450/pk) + 0.1% TEA (100mL)	µSpin-PepFrac.XX8.50
BioSPE™ PepFrac – Mini Spin	Spin columns (50/pk) + collection tubes (450/pk) + 0.1% TEA (100mL)	Spin-PepFrac.XX8.50

*Replace XX by the capacity (S or H, see p15) of your choice. e.g.: *KT-Adapt-Tips-PepFrac.H8.200.96*

BioSPE™ PepFrac for peptide fractionation

Designation	Kit composition	Reference*
BioSPE™ PepFrac – Spin 15mL	Spin columns (25/pk) + collection tubes (225/pk) + 0.1% TEA (200mL)	Spin15-PepFrac.XX8.25
BioSPE™ PepFrac – Spin 50mL	Spin columns (25/pk) + collection tubes (225/pk) + 0.1% TEA (500mL)	Spin50-PepFrac.XX8.25
BioSPE™ PepFrac – 96 plate for microelution & collection plates	Microelution 96 SPE Well Plate (1/pk) + collection plates (9/pk) + 0.1% TEA (200mL)	KT-Coll-μ96W-PepFrac.XX8.1
BioSPE™ PepFrac – 96 plate 1mL & collection plates	1mL 96 SPE Well Plate (1/pk) + collection plates (9/pk) + 0.1% TEA (500mL)	KT-Coll-96W-PepFrac.XX8.1
BioSPE™ PepFrac – Cartridges 1mL	1mL SPE cartridges (50/pk) + 0.1% TEA (100mL)	CAR1-PepFrac.XX8.50
BioSPE™ PepFrac – Cartridges 3mL	3mL SPE cartridges (50/pk) + 0.1% TEA (200mL)	CAR3-PepFrac.XX8.50
BioSPE™ PepFrac – Cartridges 6mL	6mL SPE cartridges (25/pk) + 0.1% TEA (200mL)	CAR6-PepFrac.XX8.25

*Replace XX by the capacity (S or H, see p15) of your choice. e.g.: *KT-Coll-96W-PepFrac.H8.1*



Application

Salts and interferences removal from protein mixture in top-down proteomics



Advantages

- Purification and desalting of intact proteins and large peptides (MW > 10kDa)
- Compatible with biological matrices (blood), plant samples or food matrices (milk)
- High recovery yields (> 90%)

Designation	Kit composition	Reference*
BioSPE™ PureProt – Tips 10µL Single Cell & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PureProt.SC.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PureProt.SC.10.96
BioSPE™ PureProt – Tips 10µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PureProt.XX.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PureProt.XX.10.96
BioSPE™ PureProt – Tips 200µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PureProt.XX.200.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PureProt.XX.200.96
BioSPE™ PureProt – Tips 1mL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PureProt.XX.1000.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PureProt.XX.1000.96
BioSPE™ PureProt – Micro Spin	Spin columns (50/pk) + collection tubes (100/pk)	µSpin-PureProt.XX.50
BioSPE™ PureProt – Mini Spin	Spin columns (50/pk) + collection tubes (100/pk)	Spin-PureProt.XX.50

*Replace XX by the capacity (S or H, see p15) of your choice. e.g.: µSpin-PureProt.S.96

Designation	Kit composition	Reference*
BioSPE™ PureProt – Spin 15mL	Spin columns (25/pk) + collection tubes (50/pk)	Spin15-PureProt.XX.25
BioSPE™ PureProt – Spin 50mL	Spin columns (25/pk) + collection tubes (50/pk)	Spin50-PureProt.XX.25
BioSPE™ PureProt – 96 plate for microelution & collection plates	Microelution 96 SPE Well Plate (1/pk) + collection plates (2/pk)	KT-Coll-μ96W-PureProt.XX.1
BioSPE™ PureProt – 96 plate 1mL & collection plates	1mL 96 SPE Well Plate (1/pk) + collection plates (2/pk)	KT-Coll-96W-PureProt.XX.1
BioSPE™ PureProt – Cartridges 1mL	1mL SPE cartridges (50/pk)	CAR1-PureProt.XX.50
BioSPE™ PureProt – Cartridges 3mL	3mL SPE cartridges (50/pk)	CAR3-PureProt.XX.50
BioSPE™ PureProt – Cartridges 6mL	6mL SPE cartridges (50/pk)	CAR6-PureProt.XX.50

*Replace XX by the capacity (S or H, see p15) of your choice. e.g.: Spin15-PureProt.S.50



Application

Selective enrichment of N-glycans or glycopeptides in glycoproteomics



Advantages

- Same sorbent can be used to enrich both glycopeptides or N-glycans (labelled or not)
- High recovery yields

Designation	Kit composition	Reference*
BioSPE™ GlycaClean – Tips 10µL Single Cell & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-GlycaClean.SC.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-GlycaClean.SC.10.96
BioSPE™ GlycaClean – Tips 10µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-GlycaClean.XX.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-GlycaClean.XX.10.96
BioSPE™ GlycaClean – Tips 200µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-GlycaClean.XX.200.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-GlycaClean.XX.200.96
BioSPE™ GlycaClean – Tips 1mL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-GlycaClean.XX.1000.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-GlycaClean.XX.1000.96
BioSPE™ GlycaClean – Micro Spin	Spin columns (50/pk) + collection tubes (100/pk)	µSpin-GlycaClean.XX.50
BioSPE™ GlycaClean – Mini Spin	Spin columns (50/pk) + collection tubes (100/pk)	Spin-GlycaClean.XX.50

*Replace XX by the capacity (S or H, see p15) of your choice. e.g.: Spin-GlycaClean.H.96

BioSPE™ GlycaClean for glycopeptide enrichment

Designation	Kit composition	Reference*
BioSPE™ GlycaClean – Spin 15mL	Spin columns (25/pk) + collection tubes (50/pk)	Spin15-GlycaClean.XX.25
BioSPE™ GlycaClean – Spin 50mL	Spin columns (25/pk) + collection tubes (50/pk)	Spin50-GlycaClean.XX.25
BioSPE™ GlycaClean – 96 plate for microelution & collection plates	Microelution 96 SPE Well Plate (1/pk) + collection plates (2/pk)	KT-Coll-μ96W-GlycaClean.XX.1
BioSPE™ GlycaClean – 96 plate 1mL & collection plates	1mL 96 SPE Well Plate (1/pk) + collection plates (2/pk)	KT-Coll-96W-GlycaClean.XX.1
BioSPE™ GlycaClean – Cartridges 1mL	1mL SPE cartridges (50/pk)	CAR1-GlycaClean.XX.50
BioSPE™ GlycaClean – Cartridges 3mL	3mL SPE cartridges (50/pk)	CAR3-GlycaClean.XX.50
BioSPE™ GlycaClean – Cartridges 6mL	6mL SPE cartridges (50/pk)	CAR6-GlycaClean.XX.50

*Replace XX by the capacity (S or H, see p15) of your choice. e.g.: *KT-Coll-μ96W-GlycaClean.S.1*

BioSPE™ BeadRem is used to easily and quickly filter-out magnetic beads remaining in peptide or protein samples after Single-Pot, Solid-Phase enhanced Sample Preparation (SP3) procedure, and that could plug the chromatographic column during LC-MS/MS analysis.


BioSPE™ BeadRem is available in two different versions, depending on the nature of your sample:

- **BioSPE™ BeadRem Aqua** for aqueous peptide samples containing less than 5% of organic solvent
- **BioSPE™ BeadRem Orga** for peptide samples containing more than 5% of organic solvent



Advantages

- **Efficient removal of magnetic beads** as small as **1µm diameter**
- **Fast and easy 3-step pass-through process**
- **High peptide/protein recovery**
- Possibility to **filter-out non-magnetic beads** such as **agarose beads**





Application

Removal of magnetic beads from aqueous samples containing up to 5% of organic solvent

Designation	Kit composition	Reference
BioSPE™ BeadRem Aqua – Tips 10µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (25/pk)	KT-Adapt-Tips-BeadRem-Aqua.S.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (96/pk)	KT-Adapt-Tips-BeadRem-Aqua.S.10.96
BioSPE™ BeadRem Aqua – Tips 200µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (25/pk)	KT-Adapt-Tips-BeadRem-Aqua.S.200.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (96/pk)	KT-Adapt-Tips-BeadRem-Aqua.S.200.96
BioSPE™ BeadRem Aqua – Tips 1mL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (25/pk)	KT-Adapt-Tips-BeadRem-Aqua.S.1000.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (96/pk)	KT-Adapt-Tips-BeadRem-Aqua.S.1000.96
BioSPE™ BeadRem Aqua – Micro Spin	Spin columns (50/pk) + collection tubes (50/pk)	µSpin-BeadRem-Aqua.S.50
BioSPE™ BeadRem Aqua – Mini Spin	Spin columns (50/pk) + collection tubes (50/pk)	Spin-BeadRem-Aqua.S.50
BioSPE™ BeadRem Aqua – Spin 15mL	Spin columns (25/pk) + collection tubes (25/pk)	Spin15-BeadRem-Aqua.S.25
BioSPE™ BeadRem Aqua – Spin 50mL	Spin columns (25/pk) + collection tubes (25/pk)	Spin50-BeadRem-Aqua.S.25



Application

Removal of magnetic beads from aqueous samples containing up to 5% of organic solvent

Designation	Kit composition	Reference
BioSPE™ BeadRem Aqua – 96 plate for microelution & collection plate	Microelution 96 SPE Well Plate (1/pk) + collection plate (1/pk)	KT-Coll- μ 96W-BeadRem-Aqua.S.1
BioSPE™ BeadRem Aqua – 96 plate 1mL & collection plate	1mL 96 SPE Well Plate (1/pk) + collection plate (1/pk)	KT-Coll-96W-BeadRem-Aqua.S.1
BioSPE™ BeadRem Aqua – Cartridges 1mL	1mL SPE cartridges (50/pk)	CAR1-BeadRem-Aqua.S.50
BioSPE™ BeadRem Aqua – Cartridges 3mL	3mL SPE cartridges (50/pk)	CAR3-BeadRem-Aqua.S.50
BioSPE™ BeadRem Aqua – Cartridges 6mL	6mL SPE cartridges (50/pk)	CAR6-BeadRem-Aqua.S.50



Application

Removal of magnetic beads from aqueous samples containing more than 5% organic solvent

Designation	Kit composition	Reference
BioSPE™ BeadRem Orga – Tips 10µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (25/pk)	KT-Adapt-Tips-BeadRem-Orga.S.10.25
	SPE Tips (96/pk) + adapters(96/pk) + collection tubes (96/pk)	KT-Adapt-Tips-BeadRem-Orga.S.10.96
BioSPE™ BeadRem Orga – Tips 200µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (25/pk)	KT-Adapt-Tips-BeadRem-Orga.S.200.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (96/pk)	KT-Adapt-Tips-BeadRem-Orga.S.200.96
BioSPE™ BeadRem Orga – Tips 1mL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (25/pk)	KT-Adapt-Tips-BeadRem-Orga.S.1000.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (96/pk)	KT-Adapt-Tips-BeadRem-Orga.S.1000.96
BioSPE™ BeadRem Orga – Micro Spin	Spin columns (50/pk) + collection tubes (50/pk)	µSpin-BeadRem-Orga.S.50
BioSPE™ BeadRem Orga – Mini Spin	Spin columns (50/pk) + collection tubes (50/pk)	Spin-BeadRem-Orga.S.50
BioSPE™ BeadRem Orga – Spin 15mL	Spin columns (25/pk) + collection tubes (25/pk)	Spin15-BeadRem-Orga.S.25
BioSPE™ BeadRem Orga – Spin 50mL	Spin columns (25/pk) + collection tubes (25/pk)	Spin50-BeadRem-Orga.S.25



Application

Removal of magnetic beads from aqueous samples containing more than 5% organic solvent

Designation	Kit composition	Reference
BioSPE™ BeadRem Orga – 96 plate for microelution & collection plate	Microelution 96 SPE Well Plate (1/pk) + collection plate (1/pk)	KT-Coll-μ96W-BeadRem-Orga.S.1
BioSPE™ BeadRem Orga – 96 plate 1mL & collection plate	1mL 96 SPE Well Plate (1/pk) + collection plate (1/pk)	KT-Coll-96W-BeadRem-Orga.S.1
BioSPE™ BeadRem Orga – Cartridges 1mL	1mL SPE cartridges (50/pk)	CAR1-BeadRem-Orga.S.50
BioSPE™ BeadRem Orga – Cartridges 3mL	3mL SPE cartridges (50/pk)	CAR3-BeadRem-Orga.S.50
BioSPE™ BeadRem Orga – Cartridges 6mL	6mL SPE cartridges (50/pk)	CAR6-BeadRem-Orga.S.50

Choose the perfect BioSPE™ product for your proteomic application

The following table will help you select the **BioSPE™** references that are the most adapted to your **proteomic applications**, depending on your sample volume (< 800µL) and protein/peptides amount. Please replace the **SB** in the given references by the sorbent you need for your application (e.g. KT-Adapt-Tips-PurePep.S.200.96).

		Sample volume			
		Up to 10µL	Up to 200µL	Up to 400µL	Up to 800µL
Protein/peptide amount	Up to 5µg	KT-Adapt-Tips- SB.SC.10.96	KT-Adapt-Tips- SB.S.200.96	KT-Adapt-Tips- SB.S.1000.96 KT-Coll-µ96W- SB.S.1	KT-Adapt-Tips- SB.S.1000.96
	Up to 7µg	KT-Adapt-Tips- SB.S.10.96			
	Up to 15µg	KT-Adapt-Tips- SB.H.10.96			
	Up to 25µg				
	Up to 35µg	KT-Coll-µ96W- SB.S.1	KT-Adapt-Tips- SB.H.200.96		
	Up to 50µg	KT-Coll-µ96W- SB.H.1	KT-Adapt-Tips- SB.H.200.96	KT-Adapt-Tips- SB.H.1000.96	KT-Adapt-Tips- SB.H.1000.96
	Up to 100µg		KT-Coll-µ96W- SB.H.1	KT-Coll-µ96W- SB.H.1	
	Up to 200µg	µSpin-SB.S.50	µSpin-SB.S.50	µSpin-SB.S.50	µSpin-SB.S.50
	Up to 300µg	Spin-SB.S.50	Spin-SB.S.50	Spin-SB.S.50	Spin-SB.S.50
	Up to 600µg	µSpin-SB.H.50	µSpin-SB.H.50	µSpin-SB.H.50	µSpin-SB.H.50
	Up to 900µg	Spin-SB.H.50	Spin-SB.H.50	Spin-SB.H.50	Spin-SB.H.50

Choose the perfect BioSPE™ product for your proteomic application

The following table will help you select the **BioSPE™** references that are the most adapted to your **proteomic applications**, depending on your sample volume (> 800µL) and protein/peptides amount. Please replace the **SB** in the given references by the sorbent you need for your application (e.g. KT-Adapt-Tips-PurePep.S.1000.96).

		Sample volume					
		Up to 1mL	Up to 3mL	Up to 4mL	Up to 6mL	Up to 22mL	
Protein/peptide amount	Up to 35 µg	KT-Adapt-Tips-SB.S.1000.96	CAR3-SB.S.50	Spin15-SB.S.25 CAR6-SB.S.50	CAR6-SB.S.50	Spin50-SB.S.25	
	Up to 100µg	KT-Adapt-Tips-SB.H.1000.96					
	Up to 250µg	CAR1-SB.S.50					
	Up to 500µg	KT-Coll-96W-SB.S.1					
	Up to 700µg	CAR1-SB.H.50					
	Up to 1.0mg	96W-SB.T2.1					
	Up to 1.5mg	KT-Coll-96W-SB.H.1	Spin15-SB.S.25				
	Up to 2.0mg	CAR3-SB.H.50	CAR3-SB.H.50	Spin15-SB.H.25 CAR6-SB.H.50	Spin50-SB.S.25		
	Up to 3.0mg	Spin50-SB.S.25	Spin50-SB.S.25				
	Up to 4.5mg	CAR6-SB.H.50 Spin15-SB.H.25	CAR6-SB.H.50 Spin15-SB.H.25				CAR6-SB.H.50
	Up to 9.0 mg	Spin50-SB.H.25	Spin50-SB.H.25	Spin50-SB.H.25	Spin50-SB.H.25	Spin50-SB.H.25	Spin50-SB.H.25




Application: Extraction and purification of metabolites and small molecules



Advantages

- **Desalting and concentration** of small molecules from **different matrices (urine, plasma, culture medium...)**
- **High recovery yields on a broad polarity range**
- **BioSPE™ PureOmics can run dry during SPE process**



Kit composition	Kit composition	Kit composition
BioSPE™ PureOmics – Tips 10µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PureOmics.H.10.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PureOmics.H.10.96
BioSPE™ PureOmics – Tips 200µL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PureOmics.H.200.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PureOmics.H.200.96
BioSPE™ PureOmics – Tips 1mL & Spin adapters	SPE Tips (25/pk) + adapters (25/pk) + collection tubes (50/pk)	KT-Adapt-Tips-PureOmics.H.1000.25
	SPE Tips (96/pk) + adapters (96/pk) + collection tubes (192/pk)	KT-Adapt-Tips-PureOmics.H.1000.96
BioSPE™ PureOmics – Micro Spin	Spin columns (50/pk) + collection tubes (100/pk)	µSpin-PureOmics.H.50
BioSPE™ PureOmics – Mini Spin	Spin columns (50/pk) + collection tubes (100/pk)	Spin-PureOmics.H.50
BioSPE™ PureOmics – Spin 15mL	Spin columns (25/pk) + collection tubes (50/pk)	Spin15-PureOmics.H.25
BioSPE™ PureOmics – Spin 50mL	Spin columns (25/pk) + collection tubes (50/pk)	Spin50-PureOmics.H.25

Designation	Kit composition	Reference
BioSPE™ PureOmics – 96 plate for microelution & collection plates	Microelution 96 SPE Well Plate (1/pk) + collection plates (2/pk)	KT-Coll-μ96W-PureOmics.H.1
BioSPE™ PureOmics – 96 plate 1mL & collection plates	1mL 96 SPE Well Plate (1/pk) + collection plates (2/pk)	KT-Coll-96W-PureOmics.H.1
BioSPE™ PureOmics – Cartridges 1mL	1mL SPE cartridges (50/pk)	CAR1-PureOmics.H.50
BioSPE™ PureOmics – Cartridges 3mL	3mL SPE cartridges (50/pk)	CAR3-PureOmics.H.50
BioSPE™ PureOmics – Cartridges 6mL	6mL SPE cartridges (50/pk)	CAR6-PureOmics.H.50

Choose the perfect BioSPE™ product for your metabolomic application

The following table will help you select the **BioSPE™ PureOmics** formats and references that are the most adapted to your metabolomic applications, depending on your sample volume.

		Format			
		SPE Tips	SPE Spin Columns	96 SPE Well Plates	SPE Cartridges
Sample volume	Up to 10µL	KT-Adapt-Tips-PureOmics.H.10.96	µSpin-PureOmics.H.50 Spin-PureOmics.H.50	KT-Coll-µ96W-PureOmics.H.1	CAR1-PureOmics.H.50
	Up to 200µL	KT-Adapt-Tips-PureOmics.H.200.96			
	Up to 400µL	KT-Adapt-Tips-PureOmics.H.1000.96		KT-Coll-96W-PureOmics.H.1	
	Up to 800µL				
	Up to 1mL		Spin15-PureOmics.H.25	CAR3-PureOmics.H.50	
	Up to 3mL				
	Up to 4mL	Spin50-PureOmics.H.25	CAR6-PureOmics.H.50		
	Up to 6mL				
	Up to 22mL		-		

AttractSPE® Disks

Spinnable, Automatable
High throughput HTS
Micro-SPE for microelution



Proteins / Peptides / DNA / Small molecules
Bioanalysis



Order on
www.affinisep.com

- A broad variety of sorbents
- Various formats: tips, spins, 96 well plates, SPE cartridges
- One sorbent – several capacities
- Several sorbents stacking for complex applications
(contact us for custom-made solutions)

Silica-based sorbents	Polymeric sorbents
C18 : strongly hydrophobic	SDB : polystyrene-divinylbenzene (PS-DVB) copolymer
C8 : moderately hydrophobic	HLB : SDB with hydrophilic moieties
C4 wide pore : weakly hydrophobic	RPS : Sulfonic modified SDB sorbent
Silica : highly polar	SAX : Strong anion exchange SDB
HILIC : polar	SCX : Strong cation exchange SDB


AttractSPE®Disks tools capacity

Each product from the AttractSPE®Disks range is available with **three different capacities, which are referred as T1, T2 or T3** in the product designations. This designation refers to the layer thickness of SPE membrane. For a more reproducible product, each tool contains only **ONE** layer.

Following references (p42 to p53) will have **XX indications, which have to be replaced by the capacity you require (T1, T2 or T3).**

- **T1 is A layer** having a capacity equivalent to the **thickness of ONE layer of SPE disks** (around 0.6mm)
- **T2 is A layer** having a capacity equivalent to the **thickness of TWO layers of SPE disks** (around 1.2mm)
- **T3 is A layer** having a capacity equivalent to the **thickness of THREE layers of disks** (around 1.8mm)

T1 and T2 are recommended for protein/peptide/DNA purification
T3 is recommended for small molecules extraction such as drugs in urine



Example

Capacity for 200µL tips


- T1: capacity up to 15µg
- T2: capacity up to 30µg
- T3: capacity up to 50µg

NEW! 10µL AttractSPE®Disks Tips, with a lower amount of sorbent (**T0 capacity**), for single cell-like analysis (see page 43)!



Advantages

- Several sorbents and stackings available
- Available as 10µL, 200µL and 1mL SPE Tips
- Automated packing methods with no tip-to-tip variability
- Easy and fast process by centrifugation without manual pipetting constraints



Designation	Description	200µL tips – 96/pk*	1mL tips – 96/pk*
AttractSPE® Tips C18	C18 membrane	Tips-C18.XX.200.96	Tips-C18.XX.1000.96
AttractSPE® Tips C4 Wide Pore	Wide pore C4 membrane	Tips-C4.XX.200.96	Tips-C4.XX.1000.96
AttractSPE® Tips C8	C8 membrane	Tips-C8.XX.200.96	Tips-C8.XX.1000.96
AttractSPE® Tips SDB	PS-DVB membrane	Tips-DVB.XX.200.96	Tips-DVB.XX.1000.96
AttractSPE® Tips RPS	RPS membrane	Tips-RPS-M.XX.200.96	Tips-RPS-M.XX.1000.96
AttractSPE® Tips HLB	HLB membrane	Tips-HLB.XX.200.96	Tips-HLB.XX.1000.96
AttractSPE® Tips SAX	SAX membrane	Tips-SAX.XX.200.96	Tips-SAX.XX.1000.96
AttractSPE® Tips SCX	SCX membrane	Tips-SCX.XX.200.96	Tips-SCX.XX.1000.96
AttractSPE® Tips Silica	Silica membrane	Tips-Si.XX.200.96	Tips-Si.XX.1000.96
AttractSPE® Tips HILIC	HILIC membrane	Tips-HILIC.XX.200.96	Tips-HILIC.XX.1000.96
AttractSPE® Adapters for Tips**	Spin adapters (96/pk) and collection tubes (96/pk) for tips centrifugation	Adapt-Tips.200.96	Adapt-Tips.1000.96


*Replace XX by the capacity (T1, T2, T3, see p45) of your choice. e.g.: *Tips-C18.T1.200.96*

** Tips adapters are sold only with AttractSPE® Tips – Cannot be sold separately



Advantages

- Ideal for small sample volume
- NEW! Reduced amount of sorbent (T0 capacity) for single cell-like analysis



Designation	Description	10µL tips – 96/pk *
AttractSPE® Tips C18	C18 membrane	Tips-C18.XX.10.96
AttractSPE® Tips C4 Wide Pore	Wide pore C4 membrane	Tips-C4.XX.10.96
AttractSPE® Tips C8	C8 membrane	Tips-C8.XX.10.96
AttractSPE® Tips SDB	PS-DVB membrane	Tips-DVB.XX.10.96
AttractSPE® Tips RPS	RPS membrane	Tips-RPS-M.XX.10.96
AttractSPE® Tips HLB	HLB membrane	Tips-HLB.XX.10.96
AttractSPE® Tips SAX	SAX membrane	Tips-SAX.XX.10.96
AttractSPE® Tips SCX	SCX membrane	Tips-SCX.XX.10.96
AttractSPE® Tips Silica	Silica membrane	Tips-Si.XX.10.96
AttractSPE® Tips HILIC	HILIC membrane	Tips-HILIC.XX.10.96
AttractSPE® Adapters for 10µL Tips**	Spin adapters (96/pk) and collection tubes (96/pk) for tips centrifugation	Adapt-Tips.10.96

*Replace XX by the capacity (T0, T1, T2, T3, see p45) of your choice. e.g.: *Tips-C18.T0.10.96*

** Tips adapters are sold only with AttractSPE®Tips – Cannot be sold separately

Save money while saving the planet with AttractSPE® Tips C18!

In the current environmental context, each contribution matters, that's why we now propose **960 SPE Tips in bulk packs** with economical and ecological advantages.



Advantages

- The bulk packaging represents a **reduction of more than 30% compared to the prices of 10 boxes of 96 tips**
- The new packaging **limits our impact on the environment** since it is made with **90% less plastic** and is **100% recyclable!**



DID YOU KNOW?


Our new eco packaging is **available for all AttractSPE® Tips C18**, whatever the size (10µL, 200µL, 1mL) and the capacity (T0, T1, T2, T3)

Designation	#/box	Reference*
AttractSPE® Tips C18 10µL	960 tips	Tips-C18.XX.10.960
AttractSPE® Tips C18 200µL	960 tips	Tips-C18.XX.200.960
AttractSPE® Tips C18 1mL	960 tips	Tips-C18.XX.1000.960


*Replace XX by the capacity (T0, T1, T2, T3, see p45) of your choice. e.g.: *Tips-C18.T1.200.960*
T0 capacity only available for 10µL Tips

AttractSPE® Disks Tips – StageTips

To help you select the **AttractSPE® Disks Tips** that are the most suitable for your biological applications based on your sample volume, the following table lists the maximum **loading volume** and the **elution volumes** of our different StageTips. Please note that these volumes, especially the elution volumes, are given as indicative values and depend on your applications and/or protocols.

Products		Loading volume	Elution volume
	AttractSPE® Disks Tips - 10µL	Up to 10µL	10µL
	AttractSPE® Disks Tips - 200µL	Up to 200µL	10 to 50µL
	AttractSPE® Disks Tips - 1mL	Up to 1mL	Up to 200µL


To help you select the **AttractSPE® Disks Tips** that are the most suitable for your biological applications based on your sample amount, the following table lists the **loading capacities** (maximal analytes quantities that can be retained) of our different products.

Products		Capacity* (µg)			
		T0	T1	T2	T3
	AttractSPE® Disks Tips - 10µL	3	7	15	25
	AttractSPE® Disks Tips - 200µL	-	15	30	50
	AttractSPE® Disks Tips - 1000µL	-	35	70	105

*Please note that these capacities are given as indicative values, since they are determined for C18 sorbents and are therefore expected to be higher for polymeric sorbents. T0 capacity only available for 10µL tips.

AttractSPE® Disks 96 well plates for microelution

AttractSPE® Disks 96 plate for microelution is a spinnable and automatable 96 SPE plate with wells similar to stagetips for a high throughput clean-up combined with microelution volumes.



Designation	Description	Reference – 1/pk *
AttractSPE® Disks 96 plate C18	C18 membrane	μ96W-C18.XX.1
AttractSPE® Disks 96 plate C4 Wide Pore	Wide pore C4 membrane	μ96W-C4.XX.1
AttractSPE® Disks 96 plate C8	C8 membrane	μ96W-C8.XX.1
AttractSPE® Disks 96 plate SDB	PS-DVB membrane	μ96W-DVB.XX.1
AttractSPE® Disks 96 plate RPS	RPS membrane	μ96W-RPS-M.XX.1
AttractSPE® Disks 96 plate HLB	HLB membrane	μ96W-HLB.XX.1
AttractSPE® Disks 96 plate SAX	SAX membrane	μ96W-SAX.XX.1
AttractSPE® Disks 96 plate SCX	SCX membrane	μ96W-SCX.XX.1
AttractSPE® Disks 96 plate Silica	Silica membrane	μ96W-Si.XX.1
AttractSPE® Disks 96 plate HILIC	HILIC membrane	μ96W-HILIC.XX.1

*Replace XX by the capacity (T1, T2, T3, see p48) of your choice. e.g.: μ96W-C18.T1.1


AttractSPE® Disks 96 well plates – 1mL

AttractSPE® Disks 96 Plate is a plate with 1mL SPE wells containing immobilized SPE membrane that enables a high throughput clean-up with the simultaneous preparation of 96 samples. Thanks to this small sorbent amount and a high efficiency, almost 100% of the original sample is recovered.



Advantages

- All sorbents available on catalog or on demand
- Easy handling with **automates** or **liquid handling robots**
- Compatible with a use by **centrifugation**, **positive pressure** or **aspiration**




Designation	Description	Reference – 1/pk*
AttractSPE® Disks 96 plate C18	C18 membrane	96W-C18.XX.1
AttractSPE® Disks 96 plate C4 Wide Pore	Wide pore C4 membrane	96W-C4.XX.1
AttractSPE® Disks 96 plate C8	C8 membrane	96W-C8.XX.1
AttractSPE® Disks 96 plate SDB	PS-DVB membrane	96W-DVB.XX.1
AttractSPE® Disks 96 plate RPS	RPS membrane	96W-RPS-M.XX.1
AttractSPE® Disks 96 plate HLB	HLB membrane	96W-HLB.XX.1
AttractSPE® Disks 96 plate SAX	SAX membrane	96W-SAX.XX.1
AttractSPE® Disks 96 plate SCX	SCX membrane	96W-SCX.XX.1
AttractSPE® Disks 96 plate Silica	Silica membrane	96W-Si.XX.1
AttractSPE® Disks 96 plate HILIC	HILIC membrane	96W-HILIC.XX.1


*Replace XX by the capacity (T1, T2, T3, see p48) of your choice. e.g.: 96W-C18.T1.1

AttractSPE® Disks 96 well plates 1mL and for microelution

To help you select the **AttractSPE® Disks 96 well plates** that are the most suitable for your biological applications based on your sample volume, the following table lists the maximum **loading volume** and the **elution volumes** of our different 96 well plates. Please note that these volumes, especially the elution volumes, are given as indicative values and depend on your applications and/or protocols.

Products		Loading volume	Elution volume
	AttractSPE® Disks 96 well plate for microelution	Up to 400µL	50 to 200µL
	AttractSPE® Disks 96 well plate 1mL	Up to 1mL	200 to 500µL

To help you select the **AttractSPE® Disks 96 well plates** that are the most suitable for your biological applications based on your sample amount, the following table lists the **loading capacities** (maximal analytes quantities that can be retained) of our different 96 well plates.


Products		Capacity* (µg)		
		T1	T2	T3
	AttractSPE® Disks 96 well plate for microelution	35	70	105
	AttractSPE® Disks 96 well plate 1mL	500	1000	1500

*Please note that these capacities are given as indicative values, since they are determined for C18 sorbents and are therefore expected to be higher for polymeric sorbents.

Accessories for 96 SPE well plates

Affinisep supplies accessories for the use of **AttractSPE® Disks 96 well plates**: sealing film and collection plates (0.4mL, 1mL or 2mL).

Designation	Description	Reference
Sealing film	Sealing film, 100/pk	96W-SEAL-100
96 well collection plate 0.4mL	96 round collection plate – U bottom – 0.4mL – 1/pk	96W-COLL-0.4ML-1
96 well collection plate 0.4mL	96 round collection plate – U bottom – 0.4mL – 10/pk	96W-COLL-0.4ML-10
96 well collection plate 0.4mL	96 round collection plate – U bottom – 0.4mL – 100/pk	96W-COLL-0.4ML-100
96 well collection plate 1mL	96 round collection plate – U bottom – 1mL – 1/pk	96W-COLL-1ML-1
96 well collection plate 1mL	96 round collection plate – U bottom – 1mL – 10/pk	96W-COLL-1ML-10
96 well collection plate 1mL	96 round collection plate – U bottom – 1mL – 100/pk	96W-COLL-1ML-100
96 well collection plate 2mL	96 round collection plate – U bottom – 2mL – 1/pk	96W-COLL-2ML-1
96 well collection plate 2mL	96 round collection plate – U bottom – 2mL – 10/pk	96W-COLL-2ML-10
96 well collection plate 2mL	96 round collection plate – U bottom – 2mL – 100/pk	96W-COLL-2ML-100



AttractSPE® Disks Spin Columns



Advantages

- 4 spin sizes available : micro and mini spin, 15mL and 50mL spin tubes
- Fast and easy extraction process by centrifugation
- High throughput purification
- Micro and mini spins supplied with 2mL centrifuge tubes compatible with any benchtop microcentrifuge



Micro spin Mini spin

Designation	Description	Reference micro spin – 96/pk*	Reference mini spin – 96/pk*
AttractSPE® Spin C18	C18 membrane	μSpin-C18.XX.96	Spin-C18.XX.96
AttractSPE® Spin C4 Wide Pore	Wide pore C4 membrane	μSpin-C4.XX.96	Spin-C4.XX.96
AttractSPE® Spin C8	C8 membrane	μSpin-C8.XX.96	Spin-C8.XX.96
AttractSPE® Spin SDB	PS-DVB membrane	μSpin-DVB.XX.96	Spin-DVB.XX.96
AttractSPE® Spin RPS	RPS membrane	μSpin-RPS-M.XX.96	Spin-RPS-M. XX.96
AttractSPE® Spin HLB	HLB membrane	μSpin-HLB.XX.96	Spin-HLB.XX.96
AttractSPE® Spin SAX	SAX membrane	μSpin-SAX.XX.96	Spin-SAX.XX.96
AttractSPE® Spin SCX	SCX membrane	μSpin-SCX.XX.96	Spin-SCX.XX.96
AttractSPE® Spin Silica	Silica membrane	μSpin-Si.XX.96	Spin-Si.XX.96
AttractSPE® Spin HILIC	HILIC membrane	μSpin-HILIC.XX.96	Spin-HILIC.XX.96

*Replace XX by the capacity (T1, T2, T3, see p52) of your choice. e.g.: μSpin-C18.T1.96

AttractSPE® Disks Spin Columns

AttractSPE® Spin columns of 15mL and 50mL are the perfect tools to scale up your process thanks to a high area of SPE disks.




Designation	Description	Reference 15mL spin tube – 50/µk*	Reference 50mL spin tube – 50/µk*
AttractSPE® Spin C18	C18 membrane	Spin15-C18.XX.50	Spin50-C18.XX.50
AttractSPE® Spin C4 Wide Pore	Wide pore C4 membrane	Spin15-C4.XX.50	Spin50-C4.XX.50
AttractSPE® Spin C8	C8 membrane	Spin15-C8.XX.50	Spin50-C8.XX.50
AttractSPE® Spin SDB	PS-DVB membrane	Spin15-DVB.XX.50	Spin50-DVB.XX.50
AttractSPE® Spin RPS	RPS membrane	Spin15-RPS-M.XX.50	Spin50-RPS-M.XX.50
AttractSPE® Spin HLB	HLB membrane	Spin15-HLB.XX.50	Spin50-HLB.XX.50
AttractSPE® Spin SAX	SAX membrane	Spin15-SAX.XX.50	Spin50-SAX.XX.50
AttractSPE® Spin SCX	SCX membrane	Spin15-SCX.XX.50	Spin50-SCX.XX.50
AttractSPE® Spin Silica	Silica membrane	Spin15-Si.XX.50	Spin50-Si.XX.50
AttractSPE® Spin HILIC	HILIC membrane	Spin15-HILIC.XX.50	Spin50-HILIC.XX.50


*Replace XX by the capacity (T1, T2, T3, see p52) of your choice. e.g.: Spin15-C18.T3.50

AttractSPE® Disks Spin Columns

To help you select the **AttractSPE® Disks Spin columns** that are the most suitable for your biological applications based on your sample volume, the following table lists the maximum **loading volume** and the **elution volumes** of our different Spin columns. Please note that these volumes, especially the elution volumes, are given as indicative values and depend on your applications and/or protocols.

Products		Loading volume	Elution volume
	AttractSPE® Disks Micro Spin	Up to 800µL	100 to 400µL
	AttractSPE® Disks Mini Spin	Up to 800µL	100 to 400µL
	AttractSPE® Disks Spin - 15mL	Up to 4mL	Up to 2mL
	AttractSPE® Disks Spin - 50mL	Up to 22mL	Up to 10mL

To help you select the **AttractSPE® Disks Spin columns** that are the most suitable for your biological applications based on your sample amount, the following table lists the **loading capacities** (maximal analytes quantities that can be retained) of our different Spin columns.

Products	Capacity* (µg)			
	T1	T2	T3	
	AttractSPE® Disks Micro Spin	200	400	600
	AttractSPE® Disks Mini Spin	300	600	900
	AttractSPE® Disks Spin - 15mL	1500	3000	4500
	AttractSPE® Disks Spin - 50mL	3000	6000	9000

*Please note that these capacities are given as indicative values, since they are determined for C18 sorbents and are therefore expected to be higher for polymeric sorbents.

AttractSPE® Disks Cartridges



Advantages

- 1, 3 and 6mL formats
- Large loading volume with a minimal elution volume
- High extraction capacity
- Broad range of sorbents or sorbent combinations




Designation	Reference – 1mL – 50/pk*	Reference – 3mL – 50/pk*	Reference – 6mL – 50/pk*
AttractSPE® Disks Cartridge C18	CAR1- C18.XX.50	CAR3-C18.XX.50	CAR6- C18.XX.50
AttractSPE® Disks Cartridge C4 Wide Pore	CAR1-C4.XX.50	CAR3-C4.XX.50	CAR6-C4.XX.50
AttractSPE® Disks Cartridge C8	CAR1-C8.XX.50	CAR3-C8.XX.50	CAR6-C8.XX.50
AttractSPE® Disks Cartridge SDB	CAR1- DVB.XX.50	CAR3-DVB.XX.50	CAR6- DVB.XX.50
AttractSPE® Disks Cartridge RPS	CAR1-RPS- M.XX.50	CAR3-RPS- M.XX.50	CAR6-RPS- M.XX.50
AttractSPE® Disks Cartridge HLB	CAR1- HLB.XX.50	CAR3-HLB.XX.50	CAR6- HLB.XX.50
AttractSPE® Disks Cartridge SAX	CAR1- SAX.XX.50	CAR3-SAX.XX.50	CAR6- SAX.XX.50
AttractSPE® Disks Cartridge SCX	CAR1- SCX.XX.50	CAR3-SCX.XX.50	CAR6- SCX.XX.50
AttractSPE® Disks Cartridges Silica	CAR1-Si.XX.50	CAR3-Si.XX.50	CAR6-Si.XX.50
AttractSPE® Disks Cartridges HILIC	CAR1- HILIC.XX.50	CAR3- HILIC.XX.50	CAR6- HILIC.XX.50


*Replace XX by the capacity (T1, T2, T3, see p54) of your choice. e.g.: CAR1-C18.T2.50

AttractSPE®Disks Cartridges

To help you select the **AttractSPE®Disks cartridges** that are the most suitable for your biological applications based on your sample volume, the following table lists the maximum **loading volume** and the **elution volumes** of our different cartridges. Please note that these volumes, especially the elution volumes, are given as indicative values and depend on your applications and/or protocols.

Products		Loading volume	Elution volume
	AttractSPE®Disks Cartridges - 1mL	Up to 1mL	100 to 200µL
	AttractSPE®Disks Cartridges - 3mL	Up to 3mL	200 to 400µL
	AttractSPE®Disks Cartridges - 6mL	Up to 6mL	Up to 1mL

To help you select the **AttractSPE®Disks cartridges** that are the most suitable for your biological applications based on your sample amount, the following table lists the **loading capacities** (maximal analytes quantities that can be retained) of our different cartridges.

Products		Capacity* (µg)		
		T1	T2	T3
	AttractSPE®Disks Cartridges - 1mL	250	500	750
	AttractSPE®Disks Cartridges - 3mL	700	1400	2100
	AttractSPE®Disks Cartridges - 6mL	1500	3000	4500

*Please note that these capacities are given as indicative values, since they are determined for C18 sorbents and are therefore expected to be higher for polymeric sorbents.

Frequently asked questions about

AttractSPE[®] Disks

and

BioSPE[™]



Questions about the choice of format

Q. What is the difference between mini spin and micro spin columns?

The difference between mini spin and micro spin columns is the diameter of the SPE membrane immobilized in the column (7.5mm for mini spin and 5mm for micro spin). Thus, capacities are higher for mini spin columns than for micro spin columns.

Q. What is the difference between 96 plates for microelution and 1mL 96 plates?

96 plates for microelution are plates for which the SPE wells have the shape of tips, and the volume of the wells is 400µL. 1mL 96 plates are classical plates, with SPE wells having a volume of 1mL and higher capacities than 96 well plates for microelution.

Questions concerning applications

Q. Which product can I use to purify peptides after enzymatic digestion?

For peptide purification and desalting, we recommend using BioSPE™ PurePep or BioSPE™ PurePep Broad sorbents.

Q. Which product can I use to purify large peptides and/or proteins?

For desalting large peptides and proteins, we recommend using BioSPE™ PureProt sorbent.

Q. Which product should I use for single cell like analysis?

We supply 10µL StageTips with a reduced amount of sorbent ideal for low input material and single cell-like analysis. For BioSPE™ Tips, the corresponding capacity indicated in the product reference is SC.

Q. What product is recommended for the extraction of small molecules from biological samples?

BioSPE™ PureOmics is recommended to extract metabolites and small molecules such as drugs from biological matrices (urine, plasma, blood...).

Q. I use magnetic beads for SP3 protocols, what can I use to remove them before mass spectrometry analysis of my sample?

Removal of magnetic beads or non-magnetic beads prior to MS analysis can be performed using BioSPE™ BeadRem, which is used as a specific filter to retain beads while proteins/peptides pass through the sorbent.

Questions regarding the use of our products

Q. Can I elute my sample from AttractSPE®Disks and BioSPE™ 96 well plates under vacuum?

Yes, elution from 96 well plates for microelution and 1mL 96 well plates can be achieved under vacuum, positive pressure or by centrifugation.

Q. Can I use the AttractSPE®Disks and BioSPE™ 96 well plates partially?

Yes, you can use only some wells of the AttractSPE®Disks and BioSPE™ 96 well plates for microelution and 1mL 96 well plates to perform your experiments, without consequences on the unused wells. You will just have to cover the unused wells before using the plate.

Q. Can AttractSPE®Disks and BioSPE™ Tips be used by pipetting?

To obtain optimal results, we strongly recommend using AttractSPE®Disks and BioSPE™ Tips by centrifugation or with a positive pressure manifold or with automates (refilling by the top of the tips).

Q. Are the BioSPE™ Tips supplied with adapters for centrifugation and collection tubes?

The BioSPE™ Tips are supplied as ready-to-use kits with spin adapters and collection tubes.

Q. Is there any protocol of use for AttractSPE®Disks and BioSPE™ products?

BioSPE™ products are ready-to-use kits supplied with detail instruction sheets for a given application (peptide desalting, fractionation...).

For AttractSPE®Disks products, the protocols supplied are very general to cover a wide range of applications.

Questions about Affinisep technologies

Q. What are the performances of your technology compared to other SPE products on the market?

C18 AttractSPE®Disks Tips have been compared to two brands of StageTips available on the market and showed better results for the purification of tryptic digests. The results are presented in an application note on our website ([AttractSPE®Disks Tips for an efficient protein cleanup – Comparison with competitor products](#)).

Q. What is the protein recovery ratio with the SPE membrane technology?

Thanks to their multiple advantages, our SPE membranes allow the recovery of up to 100% of the samples, with low elution volumes. A recent application note, available on our website ([Estimation of the working range on BioSPE™ PurePep SPE Tips – Standard capacity](#)), has demonstrated that, for proteins amounts as low as 10ng, the recoveries were close to 100% with RSD inferior to 3%, after desalting on 200µL BioSPE™ PurePep SEP Tips. Moreover, no loss of performance was observed for the purification of 10µg of peptides on the BioSPE™ PurePep SEP Tips standard capacity, meaning that these StageTips can be used for much higher quantities of peptides.

Questions about your specific project

Q. Is it possible to buy custom-made products that are not on Affinisep catalog?

Yes, Affinisep can provide custom products that are the most adapted to your applications. Contact our technical service at contact@affinisep.com to describe your needs and we will try our best to offer you a product that meets your criteria.

Q. I have questions regarding the use of your products, how can I contact the technical support?

You can contact our technical support team by email at tech.support@affinisep.com. You can also fill the request service on our website to describe your issue or your project.

Other products for biological applications

AttractSPE® LipRem

AttractFiltr Protein Precipitation plate

AttractFiltr membrane filter

Powder-based SPE cartridges and plates



AttractSPE® LipRem

AttractSPE® LipRem is used for the removal of phospholipids from plasma for the recovery of metabolites and small molecules.




Advantages

- Easy and fast pass-through SPE method that retains phospholipids but not small analytes
- Flow-through collection by centrifugation, positive pressure or aspiration
- Compatible with aqueous samples or lipid-based solutions
- Available as SPE cartridges, 96 well plates and reversible cartridges

Format	#/box	Reference
1mL cartridges	100	LipRem-100.S.1.20
3mL cartridges	50	LipRem-50.S.3.50
6mL cartridges	50	LipRem-50.S.6.100
96 well plate	1	LipRem-1.96W.20
Reversible cartridges, 0.7mL	25	LipRem-25.REV.1.F
	50	LipRem-50.REV.1.F

AttractFiltration Protein Precipitation plate

Removal of plasma and serum proteins by precipitation and filtration on a spinnable 96 well plate compatible with SPE automates, for the recovery of small molecules and metabolites with high recovery yields.





Designation	#/box	Reference
AttractFiltration Protein precipitation	1	PPT-1.96W
	5	PPT-5.96W



Advantages

- No leak of organic solvent used to precipitate proteins (e.g. acetonitrile)
- Filtrate collection by centrifugation, positive pressure or vacuum manifold
- Ideal for high throughput experiments
- Simple and fast 3-steps procedure



AttractFiltration membrane filter cartridges, spin tubes and 96 well plates

AttractFiltration is a filtration tool based on the use of a membrane to filtrate and remove particles with vacuum/positive pressure manifolds, SPE automates or by centrifugation before LC analysis.

Available formats are 3 and 6mL cartridges, spin tubes and 96 microfilter plates.

AttractFiltration PVDF with a PVDF membrane (hydrophobic, wide chemical compatibility, temperature resistance) for the filtration of protein digests.

AttractFiltration PTFE with a PTFE membrane (hydrophobic, wide chemical compatibility, temperature resistance) for the filtration of aggressive solutions.

AttractFiltration PES with a PES membrane (hydrophilic, low protein binding) for water filtration.

AttractFiltration Nylon with a Nylon membrane (hydrophilic, high protein, RNA & DNA binding, high surface area) for a wide range of biological preparations.

AttractFiltration CA with a Cellulose Acetate membrane (hydrophilic, low protein binding) for protein filtration.

Designation	Membrane	Pore size* (µm)	3mL Cartridges (100/pk)	6mL Cartridges (100/pk)	Spin tubes (100/pk)	96 filter plate (1/pk)
AttractFiltration PVDF	PVDF	0.45	PVDF-100.S.3.45	PVDF-100.S.6.45	PVDF-100.SPIN.45	PVDF-1.96W.45
AttractFiltration PTFE	PTFE	0.45	PTFE-100.S.3.45	PTFE-100.S.6.45	PTFE-100.SPIN.45	PTFE-1.96W.45
AttractFiltration PES	PES	0.45	PES-100.S.3.45	PES-100.S.6.45	PES-100.SPIN.45	PES-1.96W.45
AttractFiltration Nylon	Nylon	0.45	NY-100.S.3.45	NY-100.S.6.45	NY-100.SPIN.45	NY-1.96W.45
AttractFiltration CA	Cellulose acetate	0.45	CA-100.S.3.45	CA-100.S.6.45	CA-100.SPIN.45	CA-1.96W.45


*Porosity of 0.22µm also available, please contact us for more references.

Powder-based SPE cartridges and 96 well plates

AttractSPE® Polymeric-based SPE

Advantages

- **1, 3, 6, 10, 12 and 20mL cartridges** (10mg-2g of sorbent)
- **1 or 2mL 96 well plates** (10-60mg of sorbent)
- **A broad range of sorbents** (HLB, SAX, SCX, DVB, WAX, WCX)



Sorbent*	Format, sorbent amount*	#/box	40µm particles
AttractSPE® HLB	1mL cartridge, 10mg	100	HLB-100.S.1.10
	3mL cartridge, 60mg	50	HLB-50.S.3.60
	96 well plate, 10mg	1	HLB-1.96W.10
	96 well plate, 30mg	1	HLB-1.96W.30
	96 well plate, 60mg	1	HLB-1.96W.60
AttractSPE® SAX	1mL cartridge, 10mg	100	SAX-100.S.1.10
	3mL cartridge, 60mg	50	SAX-50.S.3.60
	96 well plate, 10mg	1	SAX-1.96W.10
	96 well plate, 30mg	1	SAX-1.96W.30
	96 well plate, 60mg	1	SAX-1.96W.60

*Other formats, sorbents and sorbent amounts available, please contact us.

Powder-based SPE cartridges and 96 well plates

SilactSPE Inorganic-based SPE

SilactSPE products are silica-based and alumina-based sorbents available in different formats including SPE cartridges, spin columns and 96 well plates.



Advantages

- **1, 3, 6, 10, 12mL cartridges** (50mg-2g of sorbent)
- **96 well plates** (50-100mg of sorbent)
- **A broad range of sorbents** (C18 100Å, C18 300Å, C8, C4 Wide Pore, Silica, SiSCX, SiSAX, Florisil, Phenyl, Cyano, SiWAX, SiWCX)



Sorbent*	Format, sorbent amount*	#/box	Reference
SilactSPE C4 Wide Pore	3mL cartridge, 200mg	50	C4WP-50.S.3.200
	6mL cartridge, 500mg	50	C4WP-50.S.6.500
	96 well plate, 50mg	1	C4WP-1.96W.50
SilactSPE C18 100Å	1mL cartridge, 50mg	100	C18-100.S.1.50
	3mL cartridge, 200mg	50	C18-50.S.3.200
	96 well plate, 50mg	1	C18-1.96W.50
SilactSPE C18 300Å	1mL cartridge, 50mg	100	C18WP-100.S.1.50
	3mL cartridge, 200mg	50	C18WP-50.S.3.200
	96 well plate, 50mg	1	C18WP-1.96W.50
SilactSPE C8	1mL cartridge, 50mg	100	C8-100.S.1.50
	3mL cartridge, 200mg	50	C8-50.S.3.200
	96 well plate, 100mg	1	C8-1.96W.100
SilactSPE Amine (SiWAX)	1mL cartridge, 50mg	100	NH2-100.S.1.50
	3mL cartridge, 200mg	50	NH2-50.S.3.200
	96 well plate, 100mg	1	NH2-1.96W.100

*Other formats, sorbents and sorbent amounts available, please contact us.



THE ART OF MAKING SAMPLE PREPARATION EASIER

Affinisep is

- ✓ an innovative company having **world class technology portfolio for sample preparation solution** to solve bioanalytical scientists' challenges and to impact life, environment and health with science
- ✓ **fully integrated manufacturer** of comprehensive catalog of new sample preparation kits for various analysis **such as proteomics, metabolomics and bioanalysis, food safety, water analysis, ...**

Affinisep R&D and production facilities are located in Normandy. All of our products are developed and manufactured in France.





affinisep

ORDERING INFORMATION

For any order, please choose one of the following ways:

- **On our website** : www.affinisep.com
Find all our references and benefit from a **FAST** and **SIMPLE** process!
- **By email**: contact@affinisep.com

Information about new products or new matrices:

- **Newsletters**
Stay up-to-date on AFFINISEP's latest innovations by subscribing to our newsletter, delivering new information straight to your inbox! Register on our website.